

Edexcel A Biology A-Level Core Practical 4

Investigate the effect of enzyme and substrate concentrations on the initial rates of reactions.







The rate of reaction of an enzyme-controlled reaction is influenced by different factors: the temperature, pH, concentration of the substrate, and the concentration of the enzyme. The effect of each of these can be determined by changing a single variable and measuring its effect on the rate of reaction. It is important to keep all other variables constant so that they do not influence the results. Initial rate of reaction is measured because the rate of an enzyme-controlled reaction is high, this is because enzymes act as biological catalysts, so the concentration of reactants changes rapidly. The initial rate is the only point during the reaction when concentration of reactants and products is known.

Equipment

- Milk powder solution
- Trypsin solution (1%)
- Test tubes
- Test tube holder
- Stopclock
- 5cm³ pipettes
- Goggles
- Colorimeter
- Cuvettes
- Distilled water

Method

- 1. **Dilute** stock solution of trypsin with distilled water to produce solutions with concentrations of 0.2%, 0.4%, 0.6% and 0.8%.
- 1. Make a **control** by adding 2cm³ of trypsin solution and 2cm³ of distilled water. Use this to set the colorimeter absorbance to zero.
- 2. To another cuvette, add 2cm³ of milk suspension and 2cm³ of the stock trypsin solution. Mix, place in the colorimeter and measure absorbance at 15 second intervals for 5 minutes.
- 3. Rinse the cuvette with distilled water.

Rate of reaction =
$$\frac{1}{mean \ time}$$

4. Repeat step 3 at all trypsin concentrations.





Risk Assessment

Hazard	Risk	Safety Precaution	In emergency	Risk Level
Broken glass	Cuts from sharp object	Take care when handling glass objects; keep away from edge of desk	Elevate cuts; apply pressure; do not remove glass from wound; seek medical assistance	Low
Hot liquids	Scalding	Handle with care; use tongs to remove boiling tubes from water bath; wear eye protection, keep away from edge of desk	Run burn under cold water; seek medical assistance	Low
Enzymes	Allergies	Avoid contact with skin/eyes; wear eye protection	Seek assistance	Low

Graph

• Plot a graph of rate of reaction against temperature.

Conclusion

- Milk contains a white protein called casein which, when broken down, causes the milk to turn colourless. Trypsin is a protease enzyme which hydrolyses the casein.
- As concentration of trypsin increases, the number of enzyme-substrate complexes forming also increases because enzymes and substrates are more likely to collide. This means that the rate of reaction increases up to the optimum enzyme concentration.









The rate plateaus at the point where all substrates occupy an active site.
 Increasing the enzyme concentration won't increase rate as substrate concentration is limiting the rate.

Modification

The procedure can be modified to measure the effect of **substrate concentration** on initial rate of reaction by diluting the milk suspension to produce different concentrations and **controlling concentration of trypsin**.