

## WJEC Wales Biology A Level

SP 1.6a: Scientific drawing from slides of  
root tip to show stages of mitosis

Practical notes



## Introduction

**Mitosis** is a form of **cell division** that produces two genetically **identical** daughter cells. It is important for **growth**, **repair** and **replacement** as well as **asexual reproduction**.

Root tip cells are a **good specimen** to study mitosis. They contain **meristematic stem cells** (in the apical meristem) which are continuously dividing to form **specialised** cells.

## Equipment

- Onion (*Allium* sp.) with developing roots
- 1 M hydrochloric acid
- Acetic-orcein stain
- Light microscope
- Microscope slide
- Coverslip
- Scalpel
- Mounted needle
- Fine forceps
- 2× pipettes
- Watch glass
- Paper towel
- Bunsen burner

## Risk assessment

| Hazard         | Risk          | Precaution  | Emergency  |
|----------------|---------------|---|--|
| Broken glass   | Cuts          | Keep glassware away from the edge of the desk; handle microscope slides carefully | Dispose of broken glassware carefully; elevate cuts and apply pressure; do not remove glass from cuts; seek medical assistance |
| Mounted needle | Pricking skin | Hold with pointed end downwards; keep away from the edge of the desk              | Seek medical advice  |
| Biohazard      | Contamination | Cover any cuts; wash hands after handling onion; use disinfectant                 | Seek medical advice  |



|                   |                    |  |   |
|-------------------|--------------------|--|---|
| Scalpel           | Cuts               | Direction of cut away from the body; do not attempt to change blade; keep scalpel away from the edge of the desk | Elevate cuts and apply pressure; wash minor cuts in cold water; seek medical assistance       |
| Bunsen burner     | Burns              | Act sensibly around the Bunsen burner; leave to cool before handling; wear safety goggles                        | Run burn under cold water; seek medical assistance  |
| Hydrochloric acid | Irritation to skin | Wear gloves when handling HCl  | Remove contaminated clothing; run the affected area under cold water; seek medical assistance |
|                   | Irritation to eyes | Wear safety goggles  | Flood eye(s) with tap water; seek medical assistance  |
| Acetic orcein     | Skin burns         | Wear gloves when handling acetic orcein  | Remove contaminated clothing; run the affected area under cold water; seek medical assistance |
|                   | Irritation to eyes | Wear safety goggles  | Flood eye(s) with tap water; seek medical assistance  |

## Method

### Preparing the microscope slide

- Using a scalpel, cut **1 cm** from the tip of an onion root
- Use a pipette to add **2 drops** of **1 M hydrochloric acid** into the centre of a watch glass. *HCl separates the plant cells in a process known as maceration.*
- Add **20 drops** of **acetic-orcein** stain into the watch glass centre using a second pipette. *This is used to stain the cell DNA.*
- Transfer the root tip to the centre of the watch glass, placing down **tip first**. Allow to soak for approximately **6 to 7 minutes**.



- Using the scalpel, remove **2 to 3 mm** from the root tip. Place it in the centre of a microscope slide.
- Use a pipette to add **2 to 3 drops** of acetic orcein stain onto the root tip.
- Use the Bunsen burner to gently **warm** the slide for approximately **5 seconds**.
- Break up** the tissue using a mounted needle, before applying a **cover slip**. *Lower the cover slip at an angle to prevent the formation of bubbles.*
- Squash and spread the plant tissue by applying vertical thumb pressure.
- Use a paper towel to absorb any **excess** stain on the microscope slide.

### Observing the microscope slide

- Place the microscope slide under the clips on the microscope stage.
- Turn the **lowest power objective lens** on the nose piece.
- Turn the **coarse adjustment knob** to move the stage close to the lens.
- Look down the microscope and turn the **coarse adjustment knob** to **focus** the image. Turn the **fine adjustment knob** until the best image is obtained.
- Rotate to the medium power objective lens and focus using the **fine adjustment knob**.
- Locate the apical meristem (square shaped cells with large nuclei).
- Rotate to the high power objective lens and focus using the **fine adjustment knob**.
- Observe and draw cells at **interphase, prophase, metaphase, anaphase** and **telophase**.

### Tips for high-power plan drawings

- Drawing should fill at least half of the provided space
- Only draw what you can see
- Use a **sharp pencil**
- Ensure lines are **single, complete** and **non-overlapping**
- Do **not** use shading or colour
- Create **straight lines** for labels using a ruler
- Lines should **not** intersect
- Label lines should **not** have arrow heads
- Include a scale in terms of **eyepiece units**
- Include a title and objective lens power
- Include a **magnification**



## Example

Root tip cell of an onion (*allium sp.*) in anaphase (×400)

