

# **OCR (B) Chemistry A-Level**

PAG 03: Enthalpy Determination









## 3.1 Determination of enthalpy change of neutralisation

### **Equipment**


- Polystyrene cup with lid
- Thermometer
- 250 cm<sup>3</sup> glass beaker
- 1 mol dm<sup>-3</sup> NaOH solution
- 1 mol dm<sup>-3</sup> HCl solution
- Two measuring cylinders
- Stopwatch

#### Method

- 1. Pour 25 cm<sup>3</sup> of HCl into the polystyrene cup using a measuring cylinder.
- 2. Place the polystyrene cup into the 250 cm<sup>3</sup> glass beaker.
- 3. Using a different measuring cylinder, measure 25 cm<sup>3</sup> of NaOH.
- 4. Construct a table like so:

| Time / mins      | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 |
|------------------|---|---|---|---|---|---|---|---|---|---|----|
| Temperature / °C |   |   |   |   |   |   |   |   |   |   |    |

- 5. Start the stopwatch and record the initial temperature of the HCl in the cup and it's temperature for three minutes.
- 6. On the fourth minute, add the NaOH to the cup and close the lid. Do not record the temperature.
- 7. Continue recording the temperature from 5 minutes up to 10 minutes. Stir the solution constantly.



Rice University, CC BY 4.0




#### **Calculations**

- Calculate the mean temperature of the HCl over the first three minutes recorded.
- Use the maximum temperature recorded after the NaOH was added to calculate the temperature change of the reaction.
- Use Q=mcΔT to calculate the change in energy that occurs (NaOH and HCl can be said to have a density of 1g/cm³).
- Calculate the moles of water that are produced using the chemical equation for the reaction: NaOH + HCI ⇒ H<sub>2</sub>O + NaCI.
- Use  $\Delta_{neut}H = -Q/n$  to calculate the enthalpy change of the reaction.

#### **Errors**

Heat transfer to/from the surroundings.
 Use a bomb calorimeter to minimise heat transfer.



Example of a bomb calorimeter.

Rice University, CC BY 4.0

- Inaccurate temperature measurement.
  Use a digital thermometer.
- Inaccurate determination of maximum temperature.
  Repeat the experiment with measurements at closer time intervals / attach digital thermometer to data logger.







## **Risk Assessment**

| Hazard | Risk                                                                              | Correction                          |
|--------|-----------------------------------------------------------------------------------|-------------------------------------|
| NaOH   | Can cause severe eye damage and skin burns.                                       | Wear protective gloves and goggles. |
| HCI    | Corrosive, can cause skin irritation at this concentration. Damaging to the eyes. | Wear protective gloves and goggles. |