

OCR (B) Chemistry A-Level

PAG 12: Research skills



12.1 Investigating iron tablets

Equipment list


- Weighing boat
- Balance
- Five iron tablets
- Mortar and pestle
- 100 cm³ graduated volumetric flask
- Stand and clamp
- Burette
- 100 cm³ conical flask
- 250 cm³ conical flask
- 25 cm³ pipette
- Filter funnel and paper
- 150 cm³ (0.02 mol dm⁻³) potassium manganate solution
- 150 cm³ (1 mol dm⁻³) sulfuric acid
- Distilled water

Method

1. Using the mortar and pestle grind 5 iron tablets, place them on the weighing boat and measure them accurately.
2. Add the iron tablets into the 100 cm³ conical flask and add 50 cm³ sulfuric acid.
3. Stopper the conical flask and shake the contents well until the tablet is **completely** dissolved.
4. Leave it for a while until the residue has settled.
5. Without disturbing the residue, carefully filter the solution directly into the volumetric flask.
6. Rinse the residue in the filter paper into the volumetric flask using a small volume of distilled water.
7. Add the sulfuric acid to make the solution in the flask up to the mark.
8. Mix it well to ensure all the contents are mixed. This is acidified solution of iron (II) sulfate.
9. Fill the burette with potassium manganate solution.
10. Using the pipette measure out 25 cm³ of the acidified iron (II) sulfate solution. Add this into a 250 cm³ conical flask.
11. Measure 25 cm³ sulfuric acid and mix it with the solution in the conical flask.
12. Titrate this by adding potassium manganate from the burette until the first **permanent pink** colour appears.
13. Repeat the experiment until two concordant results are obtained.
14. Record the results in a suitable format and calculate the mean titre.

A key point of this practical is to research a method yourself and then correctly cite your sources of information. The method above is just an example of the kind of practical you will be carrying out. When citing your sources make sure you include the date and time accessed, along with a full URL /if using a book you could use the Harvard referencing format.





Calculations

- Work out the moles of MnO_4^- that react.
- Find out the moles of Fe^{2+} in 25 cm^3 of the standard solution.
- Find the number of moles of Fe^{2+} in 250 cm^3 which is the moles of Fe^{2+} in five tablets.
- Work out the percentage mass of FeSO_4 in one tablet.

Errors

- Ensure all of the iron tablets are dissolved
Warming using a water bath may help.
- Finding the end point of the reaction may be tricky and hence measured incorrectly
Use a white tile under the conical flask to see the colour change more prominently.
- Some Iron tablet residue was left on the weighing boat.
Reweigh the weighing boat after adding the tablets to the conical flask.

Safety

- Sulfuric acid - may irritate skin and eyes.
- Potassium manganate - harmful if swallowed.
- Wear protective goggles throughout the experiment.

