

WJEC (Eduqas) Chemistry A-level

SP C2.1a - Preparation of a Soluble Salt by Titration

Methods and images taken from the [WJEC practical handbook](#)

This work by [PMT Education](#) is licensed under [CC BY-NC-ND 4.0](#)





SP C2.1a - Preparation of a Soluble Salt by Titration

Aim

To **prepare crystals** of sodium chloride by using **titration** followed by **evaporation**.

Apparatus and Chemicals

- 50 cm³ burette
- Funnel
- 25 cm³ pipette and filler
- 100 cm³ conical flasks
- Evaporating basin
- 0.1 mol dm⁻³ NaOH solution
- 0.1 mol dm⁻³ HCl solution
- Phenolphthalein indicator

Safety Considerations

- ★ 0.1 mol dm⁻³ NaOH solution - irritant
- ★ 0.1 mol dm⁻³ HCl solution - irritant
- ★ Phenolphthalein indicator - flammable



Method

1. Using a **pipette**, measure 25.00 cm³ of NaOH solution and pour it into a conical flask.
2. Add two drops of phenolphthalein.
3. Pour the HCl solution into a burette.
4. Record the initial volume of HCl solution.
5. Add the HCl from the burette into the conical flask a little at the time while **swirling** the conical flask.
6. When the phenolphthalein starts to turn from **pink to colourless**, add the HCl solution a drop at a time until one drop is sufficient to turn the solution colourless.
7. Record the volume of HCl solution needed.
8. Carry out the titration again using 25.00 cm³ of NaOH solution and exactly the same volume of HCl solution.
9. **Do not add the indicator** this time.
10. Gently heat the solution from the conical flask in an **evaporating basin** until its volume decreases by around a half.
11. Leave the evaporating basin to **cool** allowing **crystals** to form.

