

# WJEC (Eduqas) Chemistry A-level

## SP C2.3b - Study of an 'Iodine Clock' Reaction

Methods and images taken from the [WJEC practical handbook](#)

This work by [PMT Education](#) is licensed under [CC BY-NC-ND 4.0](#)





## SP C2.3b - Study of an 'Iodine Clock' Reaction

### Aim

To study the **kinetics** of the **oxidation of iodide ions** by hydrogen peroxide in acid solution.

### Apparatus and Chemicals

- Deionised water
- Stopwatch
- 4 x 10 cm<sup>3</sup> measuring cylinders
- 25 cm<sup>3</sup> bulb/volumetric pipette with safety filler
- 50 cm<sup>3</sup> burette and funnel
- Burette clamp and stand
- 5 x 250 cm<sup>3</sup> conical flasks
- 0.1 mol dm<sup>-3</sup> H<sub>2</sub>O<sub>2</sub> solution
- 1.0 mol dm<sup>-3</sup> H<sub>2</sub>SO<sub>4</sub> solution
- 0.1 mol dm<sup>-3</sup> KI solution
- 0.005 mol dm<sup>-3</sup> Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>
- Starch solution

### Safety Considerations

- ★ 1.0 mol dm<sup>-3</sup> H<sub>2</sub>SO<sub>4</sub> solution - irritant



### Planning

1. Decide what volumes of H<sub>2</sub>O<sub>2</sub> solution and deionised water you will mix together to get at least **5 different concentrations** of H<sub>2</sub>O<sub>2</sub>. The total volume must not exceed 5 cm<sup>3</sup>.



## Method

1. In separate conical flasks, make up the **solutions** according to the table below. Do not add the hydrogen peroxide yet.
2. Add 1 cm<sup>3</sup> starch solution to each flask and **mix** thoroughly.
3. Rapidly add the hydrogen peroxide to flask 1, starting the **stopwatch** immediately after this addition.
4. **Swirl** to mix the reaction mixture thoroughly.
5. Stop timing when the solution turns **blue-black** and record the time.
6. **Repeat** the experiment for flasks 2-5.

Flask	Volume H <sub>2</sub> SO <sub>4</sub> (cm <sup>3</sup> )	Volume Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub> (cm <sup>3</sup> )	Volume KI (cm <sup>3</sup> )	Volume H <sub>2</sub> O (cm <sup>3</sup> )	Volume H <sub>2</sub> O <sub>2</sub> (cm <sup>3</sup> )
1	10	10	25		
2	10	10	25		
3	10	10	25		
4	10	10	25		
5	10	10	25		

