

# WJEC (Wales) Chemistry A-level

# SP 1.7a - Preparation of a Soluble Salt by Titration

Methods and images taken from the WJEC practical handbook

This work by PMT Education is licensed under CC BY-NC-ND 4.0











# SP 1.7a - Preparation of a Soluble Salt by Titration

#### Aim

To prepare crystals of sodium chloride by using titration followed by evaporation.

### **Apparatus and Chemicals**

- 50 cm<sup>3</sup> burette
- Funnel
- 25 cm<sup>3</sup> pipette and filler
- 100 cm<sup>3</sup> conical flasks
- Evaporating basin
- 0.1 mol dm<sup>-3</sup> NaOH solution
- 0.1 mol dm<sup>-3</sup> HCl solution
- Phenolphthalein indicator

## **Safety Considerations**

- ★ 0.1 mol dm<sup>-3</sup> NaOH solution irritant
- ★ 0.1 mol dm<sup>-3</sup> HCl solution irritant
- ★ Phenolphthalein indicator flammable



#### Method

- 1. Using a pipette, measure 25.00 cm<sup>3</sup> of NaOH solution and pour it into a conical flask.
- 2. Add two drops of phenolphthalein.
- 3. Pour the HCl solution into a burette.
- 4. Record the initial volume of HCl solution.
- 5. Add the HCl from the burette into the conical flask a little at the time while swirling the conical flask.
- 6. When the phenolphthalein starts to turn from **pink to colourless**, add the HCl solution a drop at a time until one drop is sufficient to turn the solution colourless.
- 7. Record the volume of HCl solution needed.
- 8. Carry out the titration again using 25.00 cm<sup>3</sup> of NaOH solution and exactly the same volume of HCl solution.
- 9. Do not add the indicator this time.
- 10. Gently heat the solution from the conical flask in an **evaporating basin** until its volume decreases by around a half.
- 11. Leave the evaporating basin to **cool** allowing **crystals** to form.







