

WJEC (Wales) Chemistry A-level

SP 3.2a - Simple Redox Titration

Methods and images taken from the [WJEC practical handbook](#)

This work by [PMT Education](#) is licensed under [CC BY-NC-ND 4.0](#)



SP 3.2a - Simple Redox Titration

Aim

To determine the **relative molecular mass** of an iron(II) salt by **titration** with **standard solution** of potassium manganate(VII).

Apparatus and Chemicals

- Access to 3 decimal place mass balance (minimum 2 decimal place)
- Safety goggles
- 50 cm³ burette and funnel
- 25 cm³ pipette and filler
- 250 cm³ conical flasks
- 250 cm³ volumetric flask
- Unknown iron(II) salt
- 1.0 mol dm⁻³ H₂SO₄ solution
- Standardised KMnO₄ solution (approximately 0.02 mol dm⁻³)

Safety Considerations

- ★ 1.0 mol dm⁻³ H₂SO₄ solution - irritant
- ★ KMnO₄ solution - harmful, oxidising



Method

1. Weigh out **accurately**, about 9.8 g of the iron(II) salt provided and record the mass.
2. Make the salt up to 250 cm³ of a **standard solution** in H₂SO₄ solution.
3. Titrate 25 cm³ portions of this solution against the standardised KMnO₄ solution.
4. Use your results to calculate the **relative molecular mass** of the iron(II) salt.

The overall equation for the redox reaction can be used to help calculate the relative molecular mass of the iron(II) salt:

