

Edexcel B Biology A-Level

Core Practical 12

Investigate the rate of growth of microorganisms in liquid culture



Aseptic technique is used to **avoid contamination** of the sample from outside substances, such as microorganisms. This is important to get **reliable** and **repeatable** data.

One way of measuring bacterial growth is through **colorimetry**. This is a quick and relatively simple way of measuring growth. A disadvantage of measuring growth in this way is that it doesn't provide a direct count and also counts non-viable bacteria - to provide an indication of number of bacteria, colorimetry readings have to be related to bacterial count via a **calibration curve** produced using a **haemocytometer**.

Aseptic Technique

- Wipe down surfaces with **antibacterial cleaner** both **before and after** experiment.
- Use a **Bunsen burner** in the work space so that **convection currents** draw microbes away from the culture.
- **Flame the wire hoop** before using to transfer bacteria.
- **Flame the neck of any bottles** before use to prevent any bacteria entering the vessel (air moves out so unwanted organisms don't move in).
- Keep all vessels containing bacteria **open for the minimum amount of time**.
- **Close windows and doors** to limit air currents.

Equipment

- Bench disinfectant
- Paper towels
- Bunsen burner
- Colorimeter
- Cuvettes
- Measuring cylinders
- Yeast suspension
- Microscope
- Slides and coverslip
- Pipette
- Graph paper photocopied to acetate

Method

1. Use a cuvette filled with **glucose culture medium** to set the absorbance of the colorimeter to **zero**.
2. Fill a measuring cylinder with **yeast suspension**. Transfer it into a cuvette. Measure the absorbance.
3. Repeat steps 1 and 2 at least **five times over twelve hours**.
4. To relate colorimetry readings to a direct cell count, use a **haemocytometer**, or determine your microscope's field of view with a piece of graph paper photocopied



onto acetate. Stain the yeast suspension with methylene blue, add a drop to the slide and count the visible yeast cells.

- Calculate the volume of one drop by measuring the volume of ten drops. This information will enable you to calculate an **overall density** and therefore cell count.

Risk Assessment

Hazard	Risk	Safety Precaution	In emergency	Risk Level
Disinfectant	Flammable	Keep away from naked flame	Put out fire; seek assistance	Low
Biohazard	Contamination ; infection	Use disinfectant; wash hands with soap after dissection; do not incubate at human body temperature; do not open agar plate post incubation	Seek assistance	Low/medium (depends on likelihood of bacteria sample used to cause infection)
Naked flame	Fire hazard; burns	Keep away from flammable materials; tie up long hair	Put out fire; seek assistance; run burns under cold water immediately	Low

Graph

- Plot a graph of absorbance against time.

Conclusion

- Absorbance increases when **number of bacteria increases**. This makes the suspension more opaque, transmission decreases and absorbance increases.
- Therefore, absorbance follows the same trend as the bacterial growth curve:
 - Lag** (bacteria don't grow at maximum rate as they adapt to new environmental conditions)
 - Log** (exponential growth - bacteria grow at theoretical maximum rate)
 - Stationary** (as resources begin to be used up, bacterial growth is equal to bacterial death so number doesn't change)
 - Death** (many more bacteria are dying and so number decreases)

