

WJEC (Eduqas) Chemistry A-level

SP C1.6b - Identification of Unknown Solutions by Qualitative Analysis

Methods and images taken from the WJEC practical handbook

This work by PMT Education is licensed under CC BY-NC-ND 4.0











SP C1.6b - Identification of Unknown Solutions by Qualitative **Analysis**

Aim

To plan and carry out a method to identify six inorganic salts by the interactions between their solutions.

Apparatus and Chemicals

- 20 x test tube
- Test tube rack
- 5 x dropping pipette
- Solutions of the following salts, randomly labelled A–F
 - \circ Ba(NO₃)₂
 - \circ Pb(NO₃)₂
 - MgSO₄
 - KI
 - Na₂CO₃
 - \circ Zn(NO₃)₂

Safety Considerations

- \bigstar Ba(NO₃)₂ toxic
- \bigstar Pb(NO₃)₂ toxic



Planning

- 1. Construct and complete a table to show the expected observations when each solution interacts with the other five.
- 2. Record a summary of the observations expected for each salt (e.g. two white precipitates, one yellow precipitate and two 'no change').
- 3. Construct another, similar table but this time labelled with the letters A-F rather than the names of the salts.

Method

1. Test approximately 2 cm³ of each solution with a few drops of each of the other solutions in turn and record your observations in the second table.







