

AQA Chemistry GCSE

Required Practical 2

Neutralisation [Higher Tier]

Methods taken from the AQA Required Practical Handbook






Neutralisation

Determination of the reacting volumes of solutions of a strong acid and a strong alkali by titration

Aim

Carry out an investigation to find the **concentration** of a dilute sulfuric acid solution, using a sodium hydroxide solution of known concentration.



Equipment list

- 25 cm³ volumetric pipette
- Pipette filler
- 50 cm³ burette
- 250 cm³ conical flask
- Small funnel
- Clamp stand and clamp
- White tile
- 0.1 M sodium hydroxide solution
- Sulfuric acid
- Phenolphthalein indicator

Method


1. Use the pipette to measure 25cm³ of sodium hydroxide into the conical flask.
2. Place the conical flask on a white tile.
3. Fill the burette with sulphuric acid using a funnel.
4. Record the initial reading of acid in the burette.
- **Make sure to always take readings from the bottom of the meniscus.**
5. Add a 5 drops of indicator in this case phenolphthalein to the conical flask.
6. Slowly open the burette tap while swirling the conical flask.
7. Add acid drop-by-drop near the endpoint.
- **At this point the colour will start to change slightly.**
8. Close the burette when a colour change occurs in phenolphthalein.
- **The solution turns from pink to colourless.**
9. Record the final reading of acid in the burette and calculate the titre. This is the volume of acid used to neutralise the alkali.





- Repeat until you have concordant results.
- **These are within 0.1cm^3 of each other.**
- Present results in a table and calculate the mean titre discarding any anomalies when calculating the mean.
- Calculate the number of moles of sodium hydroxide used in the titration.
- In the balanced equation the ratio between sodium hydroxide and sulphuric acid is 2:1.
Therefore to find out the moles of sulphuric acid divide the moles of sodium hydroxide by 2.
- Use the formula [concentration = moles/volume (mean titre volume)] to work out the concentration of sulphuric acid.

Diagram



Safety Precautions

- Wear safety goggles when working with acids.
- Tie hair back.
- Report any broken glassware immediately.

