

OCR Chemistry A-Level PAG 04 - Qualitative analysis of ions

Flashcards

This work by PMT Education is licensed under CC BY-NC-ND 4.0











Describe the chemical test for halide ions











Describe the chemical test for halide ions

Add dilute nitric acid to react with carbonate ions so no Ag₂CO₃ forms (white solid).

Add silver nitrate. Precipitate forms:

- White: AgCl (soluble in dilute ammonia)
- Cream: AgBr (soluble in concentrated ammonia)
- Yellow: Agl (insoluble in ammonia)









How can SO_4^{2-} ions be identified?







How can SO_4^{2-} ions be identified?

Add aqueous barium nitrate. A white precipitate (BaSO₄) will form if sulfate ions are present. This precipitate is insoluble in excess dilute strong acids.









Describe the chemical test used to detect CO₃²⁻ and HCO₃⁻ ions











Describe the chemical test used to detect CO₃²⁻ and HCO₃⁻ ions

- Add aqueous acid
- If either ion is present bubbles of CO₂ will be released
- Test the gas by bubbling it through limewater via a delivery tube. If it is CO₂, the limewater will turn cloudy









What are the solubilities of the first 3 group 2 metals (Mg, Sr, Ba) in sulfate solution?







What are the solubilities of the first 3 group 2 metals (Mg, Sr, Ba) in sulfate solution?

Mg - soluble

Sr - insoluble

Ba - insoluble











What are the solubilities of the first 3 group 2 metals (Mg, Sr, Ba) in chromate solution?











What are the solubilities of the first 3 group 2 metals (Mg, Sr, Ba) in sulfate solution?

Mg - soluble

Sr - soluble

Ba - insoluble







